Assay for Total Cholesterol

Yoseph Cherinet (BSc, MSc; Asst. Prof)



Biochemistry Lab

Rationale

- Cholesterol is steroid alcohol synthesized in hepatocytes and is necessary for the production of bile acids, steroids, and cellular membranes.
- Total cholesterol comprises all of the cholesterol found in various lipoproteins (major component of LDLs and a minority component of HDLs and VLDLs)
- Total cholesterol has been used for aids in the detection of many conditions bound to metabolic disorders.

- Certain diseases are known causes of hypercholesterolemia
 - Cholestasis
 - Hypothyroidism
 - Diabetes
 - Nephrotic syndrome
- A high-cholesterol diet is another important factor that must be considered.
- Decreased cholesterol values found in hepatic atrophy or hypoplasia, hyperthyroidism

Total Cholesterol (TC) Measurement

Specimen : Serum/ Plasma

 Measurement of total cholesterol involves the use of three enzymes: cholesterol esterase (CE), cholesterol oxidase (CO) and peroxidase (POD) in a series of reactions in which cholesteryl esters are hydrolyzed, the 3-OH group of cholesterol is oxidized and H₂O₂ is quantified.

1. Cholesteryl ester + H₂O
$$\xrightarrow{\text{Cholesteryl esterase}}$$
 Cholesterol + Fatty acid
2. Cholesterol + O₂ $\xrightarrow{\text{Cholesterol oxidase}}$ Cholestenone + H₂O₂
3. H₂O₂ + Dye $\xrightarrow{\text{Peroxidase}}$ Color



Interferences

- Reducing substances such as ascorbic acid and bilirubin can interfere with measurements by consuming H_2O_2 .
- Lipemic samples generate turbidity of sample which leads to falsely elevated values (chylomicrons gives the appearance : Postprandial)



Materials

- Spectrophotometer
- Micropipettes
- Pipette tips (Small & Large)
- Test tubes with rack
- Cholesterol Monoreagent
- Cholesterol Standard Solution
- Test Sample
- Control Solutions

Procedure

	Blank	Standard	Sample
Monoreagent	1000ul	1000ul	1000ul
Standard		10ul	
Sample			10ul

Mix and measure the absorbance at 500 \pm 10 nm after incubating at +37 ^{0}c for 5 min or 10 min. at +25 ^{0}c .

Calculation :

C_{test =} <u>A test X C std</u> A std

Cstd= 200mg/dl

Normal values (mg/dl)	Dog	Cat	Herbivorous
	126-144	90-108	36-54