


Abs = ???



Assays for Total Protein



Yoseph Cherinet (BSc, MSc Asst. Prof)

Biochemistry II Lab

Rationale

- Every function in the living system depends on proteins.
- Various proteins constitute the bulk of the total amount of proteins present in blood.
- Total protein determination is often done to diagnose nutritional problems, liver disease, kidney disease etc


Methods for total protein determination.

- Biuret method
- Direct photometric method
- Dye-binding method
- Turbidimetric /Nephelometric
- Refractometry
- Reagent Strips

1. Biuret method

Principle

- Peptide bonds react with Cu^{2+} ions in alkaline solutions to form a violet colored product.
- The intensity of the color is proportional to the amount of protein present in the reaction. The reaction mixture is measured at 540nm.



Reaction of Cu^{2+} with peptide bonds of protein

Procedure : Sample either serum /plasma

| | Blank | Standard | Sample |
|----------------|--------------|-----------------|---------------|
| Biuret reagent | 1000ul | 1000ul | 1000ul |
| Standard | ----- | 20ul | ----- |
| Sample | ----- | ----- | 20ul |

Mix and incubate the tubes at 25° for 10 minutes and read the absorbance against reagent blank. The color is stable for an hour.

Calculation


$$C_{\text{test}} = \frac{A_{\text{test}} \times C_{\text{std}}}{A_{\text{std}}}$$

| | | | | | | | |
|-----------------------|---------|---------|---------|---------|---------|---------|---------|
| Normal values in g/dl | Dog | Cat | Cow | Horse | Sheep | Goat | Man |
| | 5.5-7.5 | 5.7-8.0 | 6.2-8.2 | 5.7-7.9 | 5.9-7.8 | 6.1-7.4 | 6.6-8.7 |

2. Direct photometric methods

Principle


- Test relies on direct photometric measurement of samples.
- Aromatic rings of tyrosine and tryptophan absorb UV light at 200-225 nm and 272-290 nm.



3. Dye-binding method

Principle


- Based on the ability of proteins to bind dyes such as amido black 10B and Coomassie Brilliant Blue.
- The color change produced when the dye binds to proteins provides a measure of total protein at 595 nm.



4. Turbidimetric and nephelometric methods

Principle:

- Protein in the sample is precipitated with addition of Sulfosalicylic acid or Trichloroacetic acid to produce turbidity. Degree of turbidity measured with Turbidometric or nephelometric methods.



5. Refractometry


- This method uses an instrument called refractometer. It measures the refractive index of solutions.
- Refractive index of water is 1.330, if a solute is added the refractive index increases linearly.
- In refractometry, it is assumed that the concentration of other compounds do not vary appreciably from serum to serum & that difference in their refractive index reflects primarily differences in protein.



6. Reagent strips



- Reagent strip testing for protein uses the principle of the protein error of indicators to produce a visible colorimetric reaction.



- Proteins carries charge at physiologic PH.
- As the protein concentration increases, the color progresses through various shades of green to blue.